

Digital impressions



iTero

Search for Centric Relations Dental Lab under the add preferred tab. If this option is unavailable call iTero directly to add Centric Relations Dental Lab as a preferred lab.

When the doctor scans any case, the file goes to Align Tech and iTero team will evaluate the scan. If the scan needs to be taken again, iTero will contact the doctor. If the scan is approved, it will be sent to the lab within 8 hours.



3 Shape Trios

Use the "add connections" feature to send request to digital@centricrelations.com.

We will get the request within 5 minutes. Once approved the doctor will be able to send scans.



Sirona

Use "add favorite lab" feature to add Centric Relations Dental Lab to your favorites.

Scans can be sent right away - simply choose Centric Relations Dental Lab from the drop down menu.



Carestream

Send request to digital@centricrelations.com through "add connections".

Files can be sent through Carestream portal, or upload the .STL file through DDS portal.